



Biotronix
Healthcare®

PROCEDURAL

Trays and Kits



CONTENTS

UROLOGY AND IRRIGATION TRAYS

Preconnected Foley Catheter Trays.....	4
Intermittent Catheter Trays.....	5
Foley Insertion Trays.....	6
Irrigation Syringe Trays.....	7

MINOR PROCEDURE TRAYS

Debridement Trays.....	8
Incision and Drainage Trays.....	9
Laceration Trays.....	10
Suture/Staple Removal Trays.....	13

RESPIRATORY CARE TRAYS

Tracheostomy Trays.....	14
Suction Catheter Kits.....	15

OR PATIENT PREP TRAYS

Dry Skin Prep Trays.....	18
Shave Prep Trays.....	20

VASCULAR PROCEDURE TRAYS

IV Start Kits.....	21
Dressing Change Trays.....	22

GENERAL HOSPITAL TRAYS

Enema Trays.....	23
------------------	----

WOMENS HEALTH TRAYS

Circumcision Trays.....	24
Pelvic Exam Trays.....	25



Urology and Irrigation Trays

Preconnected Foley Catheter Trays



Custom Procedural Trays

Biotronix Healthcare procedural trays provide solutions with procedure specific products saving you time and creating cost efficiencies. Custom procedural trays are available, our account executive will work closely with you to create the perfect tray which meets your unique requirements. Ask your Biotronix Healthcare account representative how best we can serve your facility's needs.

Tray Include:

- 2000ml Urinary Drainage Bag with reflex valve and sampling port (1)
- Poly Coated Drape 18"x 17" (1)
- PVC Powder Free Gloves, medium, with wallet (1)
- Specimen Container, 2OZ (1)
- Lubricating Jelly 5g sachet, Aplicare (1)
- PVP Swabstick, 3pcs/pk, Aplicare (1)
- Prefilled syringe with sterile water, 10cc (1)
- Paper Box (1)
- Compartmentalized Plastic Tray (1)
- Paper Caution Label (1)
- CSR Wrap 21" x 21", paper/poly/paper (1)

REF: UR1414FT

Case Quantity: 10/Case

DESCRIPTION

Latex Foley Catheter, 14/FR

REF: UR21014

Case Quantity: 10/Case

DESCRIPTION

Latex Foley Catheter, 16/FR

REF: UR1428FT

Case Quantity: 10/Case

DESCRIPTION

Latex Foley Catheter, 18/FR

Urology and Irrigation Trays

Intermittent Catheter Trays



Tray Include:

- Preconnected 1000ml Collection Bag with and without (1)
- 1000ml Graduated Tray with Tyvek Lid (1)
- Poly Coated Drape 18" x 17" (1)
- Fenestrated Drape 18"*13" paper/paper/poly, white (1)
- PVC Gloves, transparent, medium, with sleeve 1Pr (1)
- Lubricating Jelly 3g sachet, Aplicare (1)
- PVP Swabstik, 3ea/pk, Aplicare (1)
- Specimen Container 4OZ W/Lid (1)
- Patient Label (1)

REF: UR21014C
Case Quantity: 20/Case

DESCRIPTION

Intermittent Catheter Tray 14FR, PVC Catheter, Latex Free

REF: UR21014
Case Quantity: 20/Case

DESCRIPTION

Intermittent Catheter Tray 14FR, Red Rubber Catheter

REF: UR21014P
Case Quantity: 20/Case

DESCRIPTION

Intermittent Catheter Tray 14FR, PVC Catheter, Preconnected 1000ml Collection Bag

REF: UR21015
Case Quantity: 20/Case

DESCRIPTION

Intermittent Catheter Tray 15FR, Red Rubber Catheter

REF: UR21015P
Case Quantity: 20/Case

DESCRIPTION

Intermittent Catheter Tray 15FR, Red Rubber Catheter, Preconnected 1000ml Collection Bag

REF: UR21016C
Case Quantity: 20/Case

DESCRIPTION

Intermittent Catheter Tray 16FR, PVC Catheter, Latex Free

REF: UR21016P
Case Quantity: 20/Case

DESCRIPTION

Intermittent Catheter Tray 16FR, PVC Catheter, Preconnected 1000ml Collection Bag

REF: UR21016
Case Quantity: 20/Case

DESCRIPTION

Intermittent Catheter Tray 16FR, Red Rubber Catheter, Coude Tip

Urology and Irrigation Trays

Foley Insertion Trays



Case Quantity: 20/Case

REF #	DESCRIPTION
UR20122	<ul style="list-style-type: none">• 1000ml Graduated Tray with Tyvek Lid (1)• Fenestrated Drape 18" x 13", paper/paper/poly, White (1)• Poly Coated Drape 18"x 17" (1)• PVC Clear Glove, medium, with pallet (1)• Lubricating Jelly 3g sachet, Aplicare (1)• PVP Swabstick, 3ea/pk, Aplicare (1)• Prefilled Syringe with sterile water, 10cc (1)• Sterile

Urology and Irrigation Trays

Bulb Irrigation Syringe Trays



REF: URHT18

Case Quantity: 20/Case

DESCRIPTION

- Plastic Tray 800ml with Tyvek Lid (1)
- Poly Coated Drape 18"x 17" (1)
- Bulb Syringe 60cc with protective cap (1)
- Graduated Jar 500cc (1)
- Alcohol Prep Pad (1)
- Sterile

Piston Syringe Irrigation Tray

REF: URHT19

Case Quantity: 20/Case

DESCRIPTION

- Plastic Tray 800ml with Tyvek Lid (1)
- Poly Coated Drape 18"x 17" (1)
- Piston Syringe 60cc with protective cap (1)
- Graduated Jar 500cc (1)
- Alcohol Prep Pad (1)
- Sterile



Minor Procedure Trays

Debridement Trays Premium



REF: OR6201

Case Quantity: 20/Case

DESCRIPTION

- CSR, 20" x 20" (1)
- Vinyl Exam gloves, large, (2pr)
- Gauze, 3" x 3" (5)
- Instrument, forceps, Adson tissue, 4 ¾", 1 x 2 teeth (1)
- Instrument, scissors, Iris, 4 ½", straight (1)
- Ruler, 24", flexible (1)
- Cotton Tipped Applicator 5" (1)
- Surgical scalpel, #15, safety (1)
- Swabstick, PVP-I, single (1)
- Towel, blue, scrim, 20"x 19"(1)

Debridement Trays Standard

REF: OR6202

Case Quantity: 20/Case

DESCRIPTION

- CSR, 20" x 20", (1)
- Exam gloves, large (1pr)
- Gauze, 3" x 3", 12-ply (15)
- Instrument, forceps, Adson tissue, 4 ¾", 1 x 2 teeth (1)
- Instrument, forceps, Halstead Mosquito, 5", curved (1)
- Instrument, scissors, Iris, 4 ½", straight (1)
- One-compartment tray (1)
- Plastic bag (1)
- Surgical scalpel, #15, safety (1)
- Towels, blue, scrim, 20"x 19" (2)



Minor Procedure Trays

Incision and Drainage Trays



REF: OR6101
Case Quantity: 20/Case

DESCRIPTION
• Drape, poly backed, fenestrated, 18" x 13", (1)
• Gauze, 4" x 4". (4)
• Instrument, Kelly hemostat, straight, 5 1/2", (1)
• Instrument, wire forceps, 5", (1)
• Plastic bag (1)
• PVP-I prep pad (1)
• Surgical blade, #11, with plastic handle (1)
• Alcohol Prep Pad (1)
• Towels, white, absorbent 18"x 13" (1)

Incision and Drainage Trays

REF: OR6102
Case Quantity: 20/Case

DESCRIPTION
• Drape, poly backed, fenestrated, 18"x 13", (1)
• Gauze, 4" x 4", (4)
• Instrument, wire forceps, 4", (1)
• Kelly forceps, straight, 5 1/2", (1)
• Plastic bag (1)
• PVP-I prep pad (1)
• Surgical blade, #11, with plastic handle (1)
• Alcohol Prep Pad (1)
• Towels, white, absorbent 18"x 13" (1)



Incision and Drainage Trays



REF: OR6103
Case Quantity: 20/Case

DESCRIPTION
• Drape, fenestrated, 26" x 18" (1)
• Gauze, 4" x 4" (8)
• Instrument, forceps, Adson tissue, 4 3/4", 1 x 2 teeth (1)
• Instrument, forceps, Halstead Mosquito, 5", curved (1)
• Instrument, probe, 5 1/2", with eye (1)
• Instrument, scissors, operating, 5"., straight, sharp/blunt (1)
• PVP-I prep pad (2)
• Surgical blade, #11, safety, stainless steel (1)
• Towel, blue, scrim 20"x 19" (1)

Minor Procedure Trays

Laceration Tray



REF: OR4004
Case Quantity: 20/Case

DESCRIPTION


- Compartmentalized Tray with Lid
- Poly Coated Drape 18"x 17", (1)
- Fenestrated Drape 26"*16", white
- Absorbent Towel 18"*17", white
- Gauze Sponge 4"*4"*8ply (4)
- Needle Holder Floor Grade, 5"
- Tissue Forceps with Teeth Wire Grade, 5"
- Sharp/Sharp Scissors Wire Grade, 4.5"

Laceration Tray


REF: OR4005
Case Quantity: 20/Case

DESCRIPTION

- Alcohol prep pad (1)
- Drape, 4", fenestration, poly backed (1)
- Gauze, 4" x 4". (4)
- Instrument, forceps, Adson tissue, 4 3/4", straight, 1 x 2 teeth (1)
- Instrument, needle holder, Baumgartner, 5 1/2", serrated jaw (1)
- Instrument, scissors, operating, 5 1/2", straight, sharp/blunt (1)
- PVP-I prep pad (1)
- Refuse bag with tie (1)
- Towel, white (1)



Laceration Tray




REF: OR4006
Case Quantity: 20/Case

DESCRIPTION

- Alcohol prep pad (1)
- Drape, 4", fenestration, poly backed (1)
- Gauze, 4" x 4". (4)
- Instrument, forceps, Adson tissue, 4 3/4", straight, 1 x 2 teeth (1)
- Instrument, needle holder, Baumgartner, 5 1/2", serrated jaw (1)
- Instrument, scissors, operating, 5 1/2", straight, sharp/blunt (1)
- PVP-I prep pad (1)
- Refuse bag with tie (1)
- Towel, white (1)

Minor Procedure Trays

Laceration Tray



REF: OR4007

Case Quantity: 20/Case

DESCRIPTION

- CSR wrap (1)
- Drape, 2 ¾", round fenestration, white (1)
- Gauze, 2" x 2" (5)
- Gauze, 4" x 4" (10)
- Instrument, forceps, Adson tissue, 4 ¾", straight, 1 x 2 teeth (1)
- Instrument, forceps, Halstead Mosquito 5", curved, serrated jaw (1)
- Instrument, needle holder, Webster, 5", smooth jaw, (1)
- Instrument, scissors, Iris, 4 ½" in., straight, sharp/sharp (1)
- Labels, blank, white, waterproof (2)
- Medicine cups, 2 oz. (2)
- Needle, 18G x 1 ½" (1)
- Needle, 25G x 5/8" (1)
- Needle, 27G x 1 ½" (1)
- One-compartment tray (1)
- Syringe, luer lock, 10 mL (1)
- Towels, blue, scrim (2)


Laceration Tray

REF: OR4008

Case Quantity: 20/Case

DESCRIPTION

- CSR wrap (1)
- Drape, 2 ¾", round fenestration, white (1)
- Gauze, 2" x 2" in. (5)
- Gauze, 4" x 4" (10)
- Instrument, forceps, Adson tissue, 4 ¾", straight, 1 x 2 teeth (1)
- Instrument, forceps, Halstead Mosquito, 5", curved, serrated jaw (1)
- Instrument, needle holder, Baumgartner, 5 ½", serrated jaw (1)
- Instrument, scissors, Iris, 4 ½", straight, sharp/sharp (1)
- Labels, blank, white, waterproof (2)
- Medicine cups, 2 oz. (2)
- Needle, 18G x 1 ½" (1)
- Needle, 25G x 5/8" (1)
- Needle, 27G x 1 ½" (1)
- One-compartment tray (1)
- Syringe, luer lock, 10 mL (1)
- Towels, blue, scrim (2)



Laceration Tray



REF: OR4009

Case Quantity: 20/Case

DESCRIPTION

- Alcohol prep pad (1)
- Drape, 4", fenestration, poly backed (1)
- Gauze, 4" x 4" (4)
- Instrument, forceps, Adson tissue, 4 ¾", straight, 1 x 2 teeth (1)
- Instrument, forceps, Halstead Mosquito Curved, 5", curved, serrated Jaw (1)
- Instrument, needle holder, Baumgartner, 5 ½", serrated jaw (1)
- Instrument, scissors, operating, 5 ½", straight, sharp/blunt (1)
- PVP-I prep pad (1)
- Refuse bag with tie (1)
- Towel, white (1)

Minor Procedure Trays

Laceration Tray



REF: OR4010
Case Quantity: 20/Case

DESCRIPTION

- CSR wrap (1)
- Drape, 2 3/4", round fenestration, white (1)
- Gauze, 2" x 2" (5)
- Gauze, 4" x 4" (10)
- Instrument, forceps, Adson tissue, 4 3/4", straight, 1 x 2 teeth (1)
- Instrument, forceps, Halstead Mosquito, 5", curved, serrated jaw (1)
- Instrument, needle holder, Webster, 5", smooth jaw (1)
- Instrument, scissors, Iris, 4 1/2", straight, sharp/sharp (1)
- Labels, blank, white, waterproof (2)
- Medicine cups, 2 oz. (2)
- Needle, 18G x 1 1/2" (1)
- Needle, 25G x 5/8 in. (1)
- Needle, 27G x 1 1/2" (1)
- Syringe, luer lock, 10 mL (1)
- Towels, blue, cotton (2)
- Two-compartment tray (1)

Laceration Tray

REF: OR4011
Case Quantity: 20/Case

DESCRIPTION

- Cotton tip applicators, 6 in., plastic (2)
- CSR wrap (1)
- Drape, 2 3/4", round fenestration, white (1)
- Gauze, 2 x 2 in. (5)
- Gauze, 4 x 4 in. 1(0)
- Instrument, forceps, Adson tissue, 4 3/4", straight, 1 x 2 teeth (1)
- Instrument, forceps, Halstead Mosquito, 5", curved, serrated jaw (1)
- Instrument, needle holder, Webster, 5", smooth jaw (1)
- Instrument, scissors, Iris, 4 1/2", straight, sharp/sharp (1)
- Labels, blank, white, waterproof (2)
- Medicine cups, 2 oz. (2)
- Needle, 18G x 1 1/2" (1)
- Needle, 25G x 5/8" (1)
- Needle, 27G x 1 1/2" (1)
- Syringe, luer lock, 10 mL (1)
- Three-compartment tray (1)
- Towels, blue, scrim (2)



Laceration Tray



REF: OR4012
Case Quantity: 20/Case

DESCRIPTION

- Drape 27 1/2" X 19 3/4", 2" x 3" oval fenestration, blue cloth (1)
- Gauze, 4" x 4" (10)
- Instrument, forceps, Adson tissue, 4 3/4", straight, 1 x 2 teeth (1)
- Instrument, forceps, Halstead Mosquito, 5", curved, serrated jaw (1)
- Instrument, needle holder, Webster, 5", smooth jaw, large rings, satin (1)
- Instrument, scissors, Iris, 4 1/2", curved, sharp/sharp (1)
- Needle, 20G x 1 in. (1)
- Needle, 25G x 5/8 in. (1)
- Needle, 27G x 11/4 in. (1)
- Syringe, luer lock, 10 mL (1)
- Three-compartment tray (1)
- Towel, blue, cotton (1)
- CSR wrap (1)

Minor Procedure Trays

Suture Removal Trays



REF: OR4001

Case Quantity: 50/Case

DESCRIPTION

- Compartmentalized Plastic Tray with Tyvek lid (1)
- Gauze Sponge 3" x 3" (1)
- Alcohol Prep Pad (1)
- PVP Prep Pad (1)
- 4" Metal Thumb Forceps (1)
- 4.5" Metal Littauer Scissors (1)
- Sterile

Staple Removal Trays

REF: OR4002

Case Quantity: 50/Case

DESCRIPTION

- Compartmentalized Plastic Tray with Tyvek lid (1)
- Gauze Sponge 3" x 3" (1)
- Alcohol Prep Pad (1)
- PVP Prep Pad (1)
- Metal Staple Remover (1)
- Sterile



Respiratory Care Trays

Tracheostomy Care Tray – Premium



REF: RE40021

Case Quantity: 20/Case

DESCRIPTION

- Poly Coated Drape 18" x 17" (1)
- Vinyl Gloves (1pr)
- Pre-cut Tracheostomy Sponges 4" x 4" (1)
- Sponge Gauze 4" x 4" (3)
- Trach Brush 5.5" (1)
- Pipe Cleaner 5" (3)
- Cotton Tipped Applicator 5" (3)
- Twill Tape with Sponge (1)
- Hydrogen Peroxide Solution (2)
- Two Compartment Tray (1)

Tracheostomy Care Tray Plus

REF: RE40022

Case Quantity: 20/Case

DESCRIPTION

- Poly Coated Drape 18" x 17" (1)
- Vinyl Gloves (1pr)
- Pre-cut Tracheostomy Dressing 4" x 4" (1)
- Non-Woven Gauze 4" x 4" (4)
- Trach Brush 5.5" (1)
- Pipe Cleaner 5" (4)
- Cotton Tipped Applicators 3" (2)
- Twill Tape (1)
- Two Compartment Tray w/ Removable Basin



Tracheostomy Care Tray – Standard



REF: RE40023

Case Quantity: 20/Case

DESCRIPTION

- Poly Coated Drape 18" x 17" (1)
- Vinyl Gloves (1pr)
- Pre-Cut Tracheostomy Dressing 4" x 4" (1)
- Non-Woven Gauze 4" x 4" (2)
- Trach Brush 5.5" (1)
- Pipe Cleaner 5" (2)
- Cotton Tipped Applicators 3" (2)
- Pop Basin (1)
- Twill Tape (1)
- One Compartment Tray (1)

Respiratory Care Trays

Suction Catheter Kits



Tray Include:

- PVC Gloves Transparent, medium, sleeved, (1pr)
- Paper Envelope (1)
- Pop-up solution Cup (1)
- Paper Pouch (1)

REF: AN4106

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 6FR, w/graduation tubing

REF: AN4108

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 8 FR, w/graduation tubing

REF: AN4110

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 10 FR, w/graduation tubing

REF: AN4112

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 12 FR, w/graduation tubing

REF: AN4114

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 14 FR, w/graduation tubing

REF: AN4116

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 16 FR, w/graduation tubing

REF: AN4118

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 18 FR, w/graduation tubing

Respiratory Care Trays

Suction Catheter Kits



Tray Include:

- PVC Gloves Transparent, medium, sleeved, (1pr)
- Paper Envelope (1)
- Pop-up solution Cup (1)
- Paper Pouch (1)

REF: AN4119

Case Quantity: 50/Case

DESCRIPTION

- Aeroflow Suction Catheter 6FR, w/graduation tubing

REF: AN4120

Case Quantity: 50/Case

DESCRIPTION

- Aeroflow Suction Catheter 8FR, w/graduation tubing

REF: AN4121

Case Quantity: 50/Case

DESCRIPTION

- Aeroflow Suction Catheter 10 FR, w/graduation tubing

REF: AN4122

Case Quantity: 50/Case

DESCRIPTION

- Aeroflow Suction Catheter 12 FR, w/graduation tubing

REF: AN4123

Case Quantity: 50/Case

DESCRIPTION

- Aeroflow Suction Catheter 14 FR, w/graduation tubing

REF: AN4124

Case Quantity: 50/Case

DESCRIPTION

- Aeroflow Suction Catheter 16 FR, w/graduation tubing

REF: AN4125

Case Quantity: 50/Case

DESCRIPTION

- Aeroflow Suction Catheter 18 FR, w/graduation tubing

Respiratory Care Trays

Suction Catheter Kits



Tray Include:

- PVC Gloves Transparent, medium, sleeved, (1pr)
- Paper Envelope (1)
- Paper Pouch (1)
- Tray (1)

REF: AN4126

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 6FR, w/graduation tubing

REF: AN4127

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 8 FR, w/graduation tubing

REF: AN4128

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 10 FR, w/graduation tubing

REF: AN4129

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 12 FR, w/graduation tubing

REF: AN4130

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 14 FR, w/graduation tubing

REF: AN4131

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 16 FR, w/graduation tubing

REF: AN4132

Case Quantity: 50/Case

DESCRIPTION

- Chimney Suction Catheter 18 FR, w/graduation tubing

OR Patient Prep Trays

Dry Skin Scrub Tray Premium



REF: OR5003

Case Quantity: 20/Case

DESCRIPTION

- Poly Coated Drape 18" x 17" (2)
- Blotting Towel 20" x 19" (2)
- Vinyl Gloves (1pr)
- Cotton Tip Applicators 5" (6)
- Winged Sponges (6)
- Winged Sponge Sticks (3)
- Plastic Forceps (1)
- CSR Wrap (1)
- 3 Compartment Tray (1)

Dry Skin Scrub Tray Plus

REF: OR5004

Case Quantity: 20/Case

DESCRIPTION

- Poly Coated Drape 18" x 17" (2)
- Blotting Towel 20" x 19" (2)
- Vinyl Gloves, wallet (1pr)
- Cotton Tip Applicators 5" (6)
- Winged Sponges (6)
- Winged Sponge Sticks (3)
- Plastic Forceps (1)
- CSR Wrap (1)
- 3 Compartment Tray (1)



Dry Skin Scrub Tray Standard



REF: OR5005

Case Quantity: 20/Case

DESCRIPTION

- Cotton tip applicator, 6" (2)
- CSR wrap, 28" x 28" (1)
- Exam gloves, large wallet, pair (1)
- Foam small wing sponge (4)
- Foam tray (1)
- Towel, blue, scrim (2)
- Towel, white, poly lined, 17" x 18" (2)

OR Patient Prep Trays

Dry Prep Tray – Eye



REF: OR5006

Case Quantity: 20/Case

DESCRIPTION

- Towel, white, poly backed, 17" x 18" (2)
- Cotton tip applicator, 6" (6)
- CSR wrap, 28" x 28" (1)
- Emesis basin, 700cc (1)
- Exam gloves, large, wallet, pair (1)
- Gauze, non-woven, 4" x 4", 4-ply (10)
- Instrument, forceps, plastic (1)
- Medicine cup, 2 oz. (3)
- Syringe, bulb, 2oz. (1)
- Towel, blue, scrim (2)

OR Patient Prep Trays

Shave Prep Tray



REF: OR4003

Case Quantity: 20/Case

DESCRIPTION

- Poly Coated Drape 18" x 17" (1)
- Absorbent Towel 18" x 17", paper/paper/ paper (1)
- Razor (1)
- Sponge w/castile soap (1)
- Non-woven gauze 3" x 3"*6ply (3)
- Cotton Tipped Applicator, 5" (3)
- Two compartment Tray with Tyvek lid (1)
- Sterile

Vascular Procedure Trays

IV Start Kit



REF: WO4102
Case Quantity: 50/Case

DESCRIPTION

- Tourniquet 3/4" x 18", latex free (1)
- Transparent Tape 3/4" x 18" (1)
- Transparent dressing, 2 3/8" x 2 3/4". (1)
- Gauze Sponge 2" x 2" x 8ply (2)
- Alcohol Prep Pad (1)
- PVP Prep Pad (1)
- Paper Patient Label (1)
- Sterile

IV Start Kits

REF: WO4103
Case Quantity: 50/Case

DESCRIPTION

- Tourniquet 3/4" x 18", latex free (1)
- Transparent Tape 3/4" x 18" (1)
- Transparent dressing, 2 3/8" x 2 3/4". (1)
- Gauze Sponge 2" x 2" x 8ply (2)
- Alcohol Prep Pad (2)
- Paper Patient Label (1)
- Sterile



Vascular Procedure Trays

Central Line Dressing Change Tray



REF: WO5001

Case Quantity: 50/Case

DESCRIPTION

- Poly Coated Drape 18" x 13" (1)
- Gauze Sponges, 4" x 4" Large (1)
- Gauze Sponges, 2" x 2" Small (1)
- Non-Adherent Pads 2" x 2" (1)
- Pre-cut Tracheostomy Dressing 4" x 4" (1)
- Ear Loop Mask (1)
- Vinyl Gloves (1pr)
- Transparent Tape 3/4" x 18" (1)
- Transparent Dressing 4" x 4 3/4" (1)
- Paper Measuring Tape (1)
- Povidone Iodine Swabstick (3/pk)
- Alcohol Swabstick (3/pk)
- Patient Label (1)

Dressing Change Tray – Premium

REF: WO5103

Case Quantity: 20/Case

DESCRIPTION

- Poly Coated Drape, Standard 18" x 17" (2)
- Poly Coated Drape, Large 24" x 24" (1)
- Gauze Sponges, 4" x 4" Large (4)
- Gauze Sponges, 2" x 2" Small (8)
- Safety Pin (1)
- Plastic Forceps (1)
- Wire Grade Thumb Forceps (2) 5"
- Floor Grade Littauer Scissors (1pr) 4.5"
- CSR Wrap (1)
- Three Compartment Tray



Dressing Change Tray – Standard



REF: WO5104

Case Quantity: 20/Case

DESCRIPTION

- Poly Coated Drape 18" x 17" (1)
- Gauze Sponges, Large 4" x 4" (4)
- Gauze Sponges, Small 2" x 2" (8)
- Waste Bag w/Tie (1)
- Plastic Forceps (2)
- Wire Grade Thumb Forceps 4" (1)
- Wire Grade Sharp/Blunt Scissors 4.5" (1pr)
- CSR Wrap (1)
- Three Compartment Tray

General Hospital Trays

Enema Tray



REF: GE41001
Case Quantity: 48/Case

DESCRIPTION

- Underpad/Towel 18"*13", white
- Castile Soap Packet
- Enema bag with tubing 1500ml Bag
- Polybag

REF: GE41002
Case Quantity: 50/Case

DESCRIPTION

- Underpad/Towel 18"*13", white
- Castile Soap Packet
- Enema Bucket with tubing 1500ml Bucket
- Polybag

Womens Health Tray

Circumcision Tray



REF: OR6202


Case Quantity: 20/Case

DESCRIPTION

- CSR wrap (1)
- Drape, 5/8 in. fenestration, 11 5/8" x 20" (1)
- Gauze, 4" x 4", 8-ply (10)
- Instrument, forceps, Halstead Mosquito, 5", curved (2)
- Instrument, forceps, Halstead Mosquito, 5", straight (1)
- Instrument, probe, 5 1/2", with eye (1)
- Instrument, scissors, Iris, 4 1/2", straight, sharp (1)
- Instrument, scissors, operating, 4 1/2", straight, sharp/blunt (1)
- Medicine cup, 2 oz. (1)
- Petrolatum dressing, 1" x 8", Xeroform (1)
- PVP-I swabsticks, 3/pack (1)
- Safety pins (2)
- Scalpel, #10 ((1)
- Two-compartment tray (1)

Womens Health Tray

Pelvic Exam Trays



REF: WH1002

Case Quantity: 30/Case

DESCRIPTION

- Sponge forceps, plastic, 9 ½" (1)
- Vaginal speculum, plastic, medium (1)
- CSR, 24" x 24" (1)
- Foam tray (1)
- Gauze, 4" x 4", 12-ply (5)
- Lubricating jelly, 3g (3)



Biotronix **Healthcare®**

4 Robert Speck Parkway #1501,
Mississauga, ON L4Z 1S1

289-813-0159

info@biotronixhealthcare.com

www.biotronixhealthcare.com

